

# Version Changes: 1.6.3 to 1.6.4

- [b2bOpticJobData](#) → [items](#) → [item](#)
  - + element PromoID & VoucherCode added
    - [forum discussion](#)
    - [forum discussion](#)
  - + element experimental added
    - [forum discussion](#)
  - + element sessionId added
    - [forum discussion](#)
  - + element complaintManagement added
    - + element complaintManagementCode added
    - + element referenceNo added
    - + element deliveryNotelId added
    - + element invoiceId added
      - [forum discussion](#)
- + ExchangeMetaData
- + ExchangeMetaDataApplication
- + element exchangeMetaData added
  - + element application added
    - + element id added
    - + element version added
    - + element experimental added
  - + element timestamp added
  - + element hostname added
  - + element id added
  - + element name added
  - + element patientid added
  - + element experimental added
    - [forum discussion](#)
- [b2bOpticJobData](#) → [items](#) → [item](#) → [pair](#) → [patient](#)
  - + element sensitivity added
  - + element farVision added
  - + element nearVision added
  - + element previousDesign added
    - [forum discussion](#)
  - [BiometerData](#)
    - + element biometerData added
    - + element axialLength added
    - + element centralCorneaThickness added
    - + element anteriorChamberDepth added
    - + element crystallineLensThickness added
      - [forum discussion](#)
- [b2bOpticJobData](#) → [items](#) → [item](#) → duration (attribut of insurance)
  - + attribut duration added
    - [forum discussion](#)
- [b2bOpticJobData](#) → [items](#) → [item](#) → [pair](#) → [patient](#) → [gender](#)

- + Item DIVERS added
  - [forum discussion](#)
- [b2bOpticJobData](#) → [items](#) → [item](#) → [pair](#) → [lens](#) → [centration](#)
  - element monocularCentrationDistance → maxOccurs changed to 2
    - [forum discussion](#)
- [b2bOpticJobData](#) → [items](#) → [item](#) → [pair](#) → [lens](#) → [geometry](#)
  - + attribut dimension (for Inset) added
    - [forum discussion](#)
- [b2bOpticJobData](#) → [header](#) → [productcatalog](#)
  - + enumerations name lc7 added
    - [forum discussion](#)
- [b2bOpticJobData](#) → [items](#) → [item](#) → [pair](#)
  - + element inkMarkings added
    - [forum discussion](#)
- [b2bOpticJobData](#) → [items](#) → [item](#) → [pair](#) → [frame](#) → [framematerials](#)
  - + enumerations FrameMaterials ORGANICMATERIALS added
  - + enumerations FrameMaterials TITAN added
    - [forum discussion](#)
- [b2bOpticLensorder](#) → [items](#) → [item](#) → [pair](#) → [frame](#) shape → [catalog](#)
  - + Attribute “resizeToBox” added
    - [forum discussion](#)
  - + Element box to give different size with boxWidth, boxhight on right and left
    - [forum discussion](#)
- [b2bOpticJobData](#) → [items](#) → [item](#) → [pair](#) → [patient](#) → [iProfilerData](#) → [pachyeye](#)
  - + element centralCorneaThickness added
    - [forum discussion](#)

From:  
<https://wiki.b2boptic.com/> - **wiki.b2bOptic.com**

Permanent link:  
[https://wiki.b2boptic.com/en:jobdata:version010604:changelog:ver010603\\_010604](https://wiki.b2boptic.com/en:jobdata:version010604:changelog:ver010603_010604)

Last update: **2022/10/19 09:19**

