


map (IProfilerMap)

b2bOpticJobData → items → item → pair → patient → iProfilerData → result → wavefrontCornea → map



point	
occurs	1..n
x	
type	float
occurs	1
unity	µm
description	
y	
type	float
occurs	1
unity	µm
description	
z	
type	float
occurs	1
unity	µm
description	

```

<xs:complexType name="IProfilerMap">
  <xs:sequence maxOccurs="unbounded">
    <xs:element name="point">
      <xs:complexType>
        <xs:sequence>
          <xs:element name="x" type="xs:float" />
          <xs:element name="y" type="xs:float" />
          <xs:element name="z" type="xs:float" />
        </xs:sequence>
      </xs:complexType>
    </xs:element>
  </xs:sequence>
</xs:complexType>
  
```

From:
<https://wiki.b2boptic.com/> - **wiki.b2bOptic.com**

Permanent link:
<https://wiki.b2boptic.com/en:jobdata:version010604:complextypes:iprofilermap:wavefrontcornea>

Last update: **2020/05/15 11:17**

