

# Version Changes: 1.6.3 to 1.6.4

- [b2bOpticLensorder](#) → [items](#) → [item](#)
  - + element PromoID & VoucherCode added
    - [forum discussion](#)
    - [forum discussion](#)
  - + element experimental added
    - [forum discussion](#)
  - + element sessionId added
    - [forum discussion](#)
  - + element complaintManagement added
    - + element complaintManagementCode added
    - + element referenceNo added
    - + element deliveryNotelId added
    - + element invoiceId added
      - [forum discussion](#)
- [b2bOpticLensorder](#) → [items](#) → [item](#) → [pair](#) → [patient](#)
  - + element sensitivity added
  - + element farVision added
  - + element nearVision added
  - + element previousDesign added
    - [forum discussion](#)
  - [BiometerData](#)
  - + element biometerData added
  - + element axialLength added
  - + element centralCorneaThickness added
  - + element anteriorChamberDepth added
  - + element crystallineLensThickness added
    - [forum discussion](#)
- [b2bOpticLensorder](#) → [items](#) → [item](#) → duration (attribut of insurance)
  - + attribut duration added
    - [forum discussion](#)
- [b2bOpticLensorder](#) → [items](#) → [item](#) → [pair](#) → [patientgender](#)
  - + Item DIVERS added
    - [forum discussion](#)
- [b2bOpticLensorder](#) → [items](#) → [item](#) → [pair](#) → [lens](#) → [centration](#)
  - element monocularCentrationDistance → maxOccurs changed to 2
    - [forum discussion](#)
- [b2bOpticLensorder](#) → [items](#) → [item](#) → [pair](#) → [lens](#) → [geometry](#)
  - attribut dimension (for Inset) added
    - [insetdimension](#)
    - [forum discussion](#)
- [b2bOpticLensorder](#) → [header](#) → [productCatalog](#)
  - Item lc7 added
    - [forum discussion](#)

- [b2bOpticLensorder](#) → [items](#) → [item](#) → [pair](#)
  - + element `inkMarkings` added
    - [forum discussion](#)
- [b2bOpticLensorder](#) → [items](#) → [item](#) → [pair](#) → [frame](#) → [framematerials](#)
  - + enumerations `FrameMaterials ORGANICMATERIALS` added
  - + enumerations `FrameMaterials TITAN` added
    - [forum discussion](#)
- [b2bOpticLensorder](#) → [items](#) → [item](#) → [pair](#) → [frame](#) shape → catalog
  - + Attribute `“resizeToBox”` added
    - [forum discussion](#)
  - + Element `box` to give different size with `boxWidth`, `boxheight` on right and left
    - [forum discussion](#)
- [b2bOpticLensorder](#) → [items](#) → [item](#) → [pair](#) → [patient](#) → [iProfilerData](#) → [result](#) → [pachyeye](#)
  - + element `centralCorneaThickness` added
    - [forum discussion](#)

From:

<https://wiki.b2boptic.com/> - **wiki.b2bOptic.com**

Permanent link:

[https://wiki.b2boptic.com/en:lensorder:changelog:ver010603\\_010604](https://wiki.b2boptic.com/en:lensorder:changelog:ver010603_010604)

Last update: **2024/07/09 04:39**

