


# PachyEye (IProfilerResult)

b2boptic → items → item → pair → patient → iProfilerData → result → pachyeye



| distCorneaLens              |                            |
|-----------------------------|----------------------------|
| type                        | float                      |
| unity                       | dpt                        |
| occurs                      | 1                          |
| description                 | distance Cornea Lens       |
| New! centralCorneaThickness |                            |
| type                        | float                      |
| unity                       | mm                         |
| occurs                      | 0..1                       |
| description                 | Cornea Thickness in middle |

```
<xs:complexType name="PachyEye">
  <xs:sequence>
    <xs:element name="distCorneaLens" type="xs:float"/>
    <xs:element name="centralCorneaThickness" type="xs:float"
minOccurs="0"/>
  </xs:sequence>
</xs:complexType>
```

From: <https://wiki.b2boptic.com/> - **wiki.b2bOptic.com**

Permanent link: <https://wiki.b2boptic.com/en:lensorder:version010604:complextypes:pachyeye>

Last update: **2021/06/14 13:03**

