


# map (IProfilerMap)

b2boptic → items → item → pair → patient → iProfilerData → result → wavefrontEye → map



| point       |       |
|-------------|-------|
| occurs      | 1..n  |
| x           |       |
| type        | float |
| occurs      | 1     |
| unity       | µm    |
| description |       |
| y           |       |
| type        | float |
| occurs      | 1     |
| unity       | µm    |
| description |       |
| z           |       |
| type        | float |
| occurs      | 1     |
| unity       | µm    |
| description |       |

```
<xs:complexType name="IProfilerMap">
  <xs:sequence maxOccurs="unbounded">
    <xs:element name="point">
      <xs:complexType>
        <xs:sequence>
          <xs:element name="x" type="xs:float" />
          <xs:element name="y" type="xs:float" />
          <xs:element name="z" type="xs:float" />
        </xs:sequence>
      </xs:complexType>
    </xs:element>
  </xs:sequence>
</xs:complexType>
```

From:

<http://wiki.b2boptic.com/> - **wiki.b2bOptic.com**

Permanent link:

<http://wiki.b2boptic.com/en:lensorder:version010605:complextypes:iprofilermmap:wavefronteye>

Last update: **2024/07/09 04:10**

