


# PachyEye (IProfilerResult)

b2boptic → items → item → pair → patient → iProfilerData → result → pachyeye



| distCorneaLens                     |                            |
|------------------------------------|----------------------------|
| type                               | float                      |
| unity                              | dpt                        |
| occurs                             | 1                          |
| description                        | distance Cornea Lens       |
| <b>New!</b> centralCorneaThickness |                            |
| type                               | float                      |
| unity                              | mm                         |
| occurs                             | 0..1                       |
| description                        | Cornea Thickness in middle |

```
<xs:complexType name="PachyEye">
  <xs:sequence>
    <xs:element name="distCorneaLens" type="xs:float"/>
    <xs:element name="centralCorneaThickness" type="xs:float"
minOccurs="0"/>
  </xs:sequence>
</xs:complexType>
```

From:

<http://wiki.b2boptic.com/> - **wiki.b2bOptic.com**

Permanent link:

<http://wiki.b2boptic.com/en:lensorder:version010605:complextypes:pachyeye>

Last update:

**2024/07/09 04:10**